# CinemaVision

### **New Horizons in Patient Comfort**

Transforming the patient's MRI experience, the CinemaVision is a versatile and immersive entertainment system. It provides multiple Audio/Video options from Blu-Ray, DVD, CD, and AUX inputs, while offering two-way communication and dramatically reducing MRI gradient noise. The Visor/Headset combination fits completely within the head coil and operates seamlessly inside the magnet bore with no detrimental effects on signal/noise.



Eliminate appointment cancellations and sedation as patients enjoy a stunning sight and sound experience. With LaserLink technology, installation is just a matter of minutes. Patient set-up is quick and simple. Call today to learn more!



#### Virtual Command Center A



The CinemaVision Controller manages entertainment choices for both the patient and technologist, enables two-way communication with patients, and more.

#### Ultra-Sleek and Comfortable

CinemaVision's integrated audio-video headset, tranducer and visor deliver top-quality entertainment and a virtual MRI comfort zone.

#### **Handsome High-Tech Remote**



The CinemaVision remote control puts the functionality of the entire CinemaVision system at your fingertips!



Resonance Technology, Inc.

818.882.1997 • 800.428.MRTV • www.mrivideo.com

This product is under one or more of the following US Patent: [9,787,750], [5,877,732], [5,412,419], [5,627,902], [5,432,544], [9,571,820], [9,829,547], [9,979,938].

CinemaVision

Setting New Standards in Patient Relaxation

## Making MRI Virtually Entertaining

Delivering unparalleled virtual reality, CinemaVision transports patients to a breathtaking world of immersive entertainment.

Combining vivid video and symphonic-quality audio, Resonance Technology offers the ultimate in patient comfort. CinemaVision comes with a lightweight, adjustable headset which can activate two way communication paired with a Video Visor.

Treat your patients to the highest quality multi-media experience and enjoy its top-of-the-line performance!



#### **CinemaVision Features Include:**

- Widescreen video display with 16:9 aspect ratio or standard television 4:3
- Total resolution of 240,000 pixels on each display
- 45 degrees horizontal field-of-view
- Video Visor that produces the equivalent of a dramatic 62-inch screen viewed from 5.5 feet
- Proprietary MRI LaserLink technology that transmits audio-video signals from the console to the magnet room
- Digital audio system with 30dB noise-attenuating headset with 500Hz - 10KHz frequency response
- Compatible with 3T-rated field strength magnets, ensuring compatibility with all state-of-the-art MRIs



Resonance Technology, Inc.

818.882.1997 • 800.428.MRTV • www.mrivideo.com

This product is under one or more of the following US Patent: [9,787,750], [5,877,732], [5,412,419], [5,627,902], [5,432,544], [9,571,820], [9,829,547], [9,979,938].