

HD MR Projector

Setting New Standards in Patient Comfort



Resonance Technology, Inc.'s high resolution projector for the MRI is finally available. This small compact MR Projector is easy to use. The Hi-Def projection image distracts the patient during the MRI scan procedure. Turning any MRI magnet into a personalized theater for the patient to relax and enjoy the vibrant video.

MR Projector Specifications:

- » Screen can be placed inside MRI bore
- » Resolution: 1920 x 1080 (1080p)
- » Image size: 5" - 120" (diagonally)
- » Image Flip and Mirror functions
- » Fits in any type of MRI suite
- » Adjustable projector stand
- » Pixel Refresh rate: 2ms
- » Fiber optic connection
- » Adjustable brightness
- » 500 Lumens



MR Projector with Lens ▲



Resonance Technology, Inc.

818.882.1997 • 800.428.MRTV • www.mrvideo.com

This product is under one or more of the following US Patent: [9,787,750], [5,877,732], [5,412,419], [5,627,902], [5,432,544], [9,571,820], [9,829,547], [9,979,938].